



The logo for OPNTEC, featuring the word "OPNTEC" in white, bold, sans-serif capital letters on an orange rectangular background.


PSLab Sensor Boxes

Rebekka Roßberg, Mario Behling

The PSLab.io logo, consisting of the text "PSLab.io" in white, bold, sans-serif font on a red rectangular background.

Let's make science and knowledge accessible for everyone! From entry level to advanced and industry, create your own scientific sensor boxes and experiments easily.





1. Measuring Light Intensity


- ❖ Measurement in Lux
- ❖ Using the **GL5528** photo-resistor
- ❖ Embedded within the **KY-018**'s voltage divider circuit
- ❖ Value on exponential scale

PSLab.io

Lux-Meter

#01

Light Intensity Sensor: GL5528




2. Measuring Temperature

- ❖ Using the **LM35** sensor
- ❖ Range: -55 °C to 150 °C
- ❖ Accuracy: +/- 1 °C

PSLab.io

Thermometer #01

Temperature Sensor: LM35



3. Measuring CO₂


- ❖ Measurement in particles per million (ppm)
- ❖ Using the **CCS811** sensor for CO₂e quantification:
 - Accurate measurement of green house gases
- ❖ Automatic offset of temperature and humidity
- ❖ Lower limit of 400ppm CO₂e
- ❖ Embedded on the **CJMCU-8118** circuit board
- ❖ With custom I²C driver

PSLab.io

CO₂ Sensor


#01

CO₂ Sensor: CCS811



4. Measuring Oxygen

- ❖ Measures the percentage of oxygen in the surrounding air
- ❖ Using the **A0-03** sensor
- ❖ Embedded in a custom circuit board for signal amplification




PSLab.io

O₂ Sensor

#01

Oxygen Sensor: A0-03



Stealth Design for Outdoor Usage




Basic Usage

- ❖ Plug Power Source into the PSLab Sensor Box
 - **Measurement starts**
 - Wi-Fi Hotspot for file sharing is opened
- ❖ To fetch the measurement data:
 - Connect PC/ phone to the corresponding **Wi-Fi Hotspot**
 - Access **file server**
 - Fetch **CSV** file
- ❖ Push button to trigger the safe shutdown

Connecting

- ❖ Detailed **manuals** to initialize a connection to the file server can be found in the [GitHub Repository](#)
- ❖ Currently Windows, Mac, Linux, Android, and iPhone are supported there



General Advice

- ❖ Wait a minute or two after starting the PSLab sensor box
- ❖ Ensure an **active Wi-Fi connection**
- ❖ Access the file server “anonymously”
- ❖ Choose unique connection names

Linux

- ❖ Open the file explorer
- ❖ The new device “PSLab” should automatically appear in the network section (else, refresh the page)
- ❖ Connect to this device

Windows

- ❖ Open the file explorer
- ❖ Right click on "This PC" and select "Add Network Address"
- ❖ Choose to add a user defined network address and enter "\\10.42.0.1\data" as this address
- ❖ Assign a name to this new connection

Mac

- ❖ Hit Command+K and look for "Connect to Server"
- ❖ Enter "smb://10.42.0.1/data" as the server address and click "Connect"
- ❖ Confirm the action by entering your login details


iPhone

- ❖ Open the "Files" app
- ❖ Select "Browse" > "Connect to Server"
- ❖ Enter "smb://10.42.0.1/data" as the network address and click "Connect"

Android

- ❖ Network device support depends on the specific phone type
- ❖ In case they are not supported, there are several external applications for this use case available in the Play Store or on F-Droid


Accessing the Data



The screenshot shows a text editor window with the file '18-05-2023(2).csv'. The file contains 15 rows of data, each representing a timestamp and sensor readings. The status bar at the bottom indicates 'Ln 1, Col 1', '100%', 'Windows (CRLF)', and 'UTF-8'.

18-05-2023	19:55:26	2.2650298683403194	lux	light
18-05-2023	19:55:27	2.2490474692128526	lux	light
18-05-2023	19:55:28	2.23309875282076	lux	light
18-05-2023	19:55:29	2.282026028134832	lux	light
18-05-2023	19:55:30	2.3159017021941731	lux	light
18-05-2023	19:55:31	2.820810481415256	lux	light
18-05-2023	19:55:32	3.6229173103752304	lux	light
18-05-2023	19:55:33	5.1852017303712943	lux	light
18-05-2023	19:55:34	5.294026907891021	lux	light
18-05-2023	19:55:35	4.9241904847927492	lux	light
18-05-2023	19:55:36	5.02819347529197432	lux	light
18-05-2023	19:55:37	5.29520271504272934	lux	light
18-05-2023	19:55:38	4.5018318493201485	lux	light
18-05-2023	19:55:39	3.293520752959271	lux	light
18-05-2023	19:55:40	2.7320520525963549	lux	light
18-05-2023	19:55:41	2.83513549673874621	lux	light

Up Next: Advanced Projects for Students



PSLab.io


- ❖ Accessing the Raspberry Pi operating system
- ❖ Manually running the sensor scripts
- ❖ Adjusting the **measurement interval**
- ❖ Adjusting the **flushing interval**
- ❖ Accessing the continuous measurement data stream via an **API**
- ❖ Fetching Software Updates from GIT
- ❖ Changing the **sensor type**
- ❖ ... and many more

Experimenting with Python

```
MEASURING_INTERVAL = 1 # in seconds
```

```
experiment_options = { # item_name : [function_name, unit]  
    "co2"      : [measure_co2, "ppm"],  
    "oxygen"   : [measure_oxygen, "%"],  
    "light"    : [measure_light_intensity, "  
    "temp"     : [measure_temperature, "°C"]  
}
```

```
In [41]: ax = df['T_ambient'].plot(figsize=(12, 6), grid=True, marker='^')  
ax = df['T_ambient'].rolling(window=12, center=True).mean().plot(ax=ax, grid=True, lw=2, c='r', marker='o')  
ax.set_title('Ambient temperature')  
ax.legend(['Ambient', 'Rolling window ambient']);
```



Thank you!

OPNTEC



<https://github.com/fossasia/pslab-scripts>

Twitter, Github, FB, Linkedin: @pslabio @opntec

PSLab.io